

## Product Datasheet

### Furin Inhibitor II (orb654641)

<b>Description</b>	Inhibitor of furin (Ki values are 0.106, 0.58 and 13.2 $\mu$ M for furin, PACE4 and PC1 respectively).
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>36</sub> H <sub>75</sub> N <sub>25</sub> O <sub>6</sub>
<b>MW</b>	954.16
<b>CAS Number</b>	[673202-67-0]
<b>SMILES</b>	[H]N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N[C@H](CCCNC(N)=N)C(=O)N
<b>Expiration Date</b>	12 months from date of receipt.