

Product Datasheet

Pexmetinib(ARRY-614) (orb654669)

| | |
|------------------------|---|
| Description | Pexmetinib is a potent Tie-2 and p38 MAPK dual inhibitor, with IC50s of 1 nM, 35 nM and 26 nM for Ti... |
| Note | For research use only |
| Formula | $C_{31}H_{33}FN_6O_3$ |
| MW | 556.25 |
| CAS Number | [945614-12-0] |
| SMILES | <chem>O=C(NCC1=CC(F)=CC=C1OC2=CC3=C(N(CCO)N=C3)C=C2)NC4=CC(C)(C)(C)=NN4C5=CC=C(C)C=C5</chem> |
| Expiration Date | 12 months from date of receipt. |