

Product Datasheet

PIK-75 free base (orb654698)

Description	PIK-75 is a selective cell-permeable inhibitor of PI 3-kinase p110 α (IC50 values are 5.8, 76 nM and...
Note	For research use only
Formula	C ₁₆ H ₁₄ BrN ₅ O ₄ S
MW	452.28
CAS Number	[372196-67-3]
SMILES	CC1=C(C=C(C=C1)[N+](=O)[O-])S(=O)(=O)N(C)/N=C\C2=CN=C3N2C=C(C=C3)Br
Expiration Date	12 months from date of receipt.