

Product Datasheet

PKC-theta inhibitor (orb654699)

Description	PKC-theta inhibitor is a selective PKC- θ inhibitor, with an IC50 of 12 nM.
Note	For research use only
Formula	$C_{20}H_{25}F_3N_6O_3$
MW	454.45
CAS Number	[736048-65-0]
SMILES	<chem>FC(F)(F)OC1=CC=CC=C1CNC2=NC=C([N+])([O-])=O)C(NCC3CCC(CN)CC3)=N2</chem>
Expiration Date	12 months from date of receipt.