

Product Datasheet

PXS-5153A (orb654705)

Description	PXS-5153A (PXS5153A) is a potent, selective, fast-acting, dual LOXL2/LOXL3 inhibitor with IC50 of < ...
Note	For research use only
Formula	$C_{20}H_{25}Cl_2FN_4O_2S$
MW	475.404
CAS Number	[2125956-82-1]
SMILES	<chem>C(S(N(C)C)(=O)=O)1=CC=CC(C2C3=NC=CC=C3N(C/C(/F)=C/CN)C=2C)=C1.[H]Cl.[H]Cl</chem>
Expiration Date	12 months from date of receipt.