

Product Datasheet

ML367 (orb654712)

Description	ML367 is a potent inhibitor of ATPase family AAA domain-containing protein 5 (ATAD5) stabilization, ...
Note	For research use only
Formula	$C_{19}H_{12}F_2N_4$
MW	334.32
CAS Number	[381168-77-0]
SMILES	<chem>FC1=CC(NC2=NC(C3=CC=NC=C3)=NC4=CC=CC=C42)=CC=C1F</chem>
Expiration Date	12 months from date of receipt.