

Product Datasheet

Israpafant (orb654754)

Description	Israpafant (Y24180) is a specific antagonist of PAF receptor that inhibits PAF-induced rabbit platelet...
Note	For research use only
Formula	$C_{28}H_{29}ClN_4S$
MW	489.074
CAS Number	[117279-73-9]
SMILES	<chem>CC1C2=NN=C(N2C3=C(C=C(S3)CCC4=CC=C(C=C4)CC(C)C(=N1)C5=CC=CC=C5Cl)C</chem>
Expiration Date	12 months from date of receipt.