

## Product Datasheet

# Human Novel Coronavirus Nucleoprotein/Spike Glycoprotein S1 (SARS-CoV-2 N/S1) IgG Antibody ELISA Kit (orb654833)

<b>Description</b>	Human Receptor for advanced glycation end products, Novel Coronavirus Nucleoprotein/Spike Glycoprote...
<b>Reactivity</b>	Human
<b>Range</b>	Request Information
<b>Target</b>	SARS-CoV-2 N/S1 Ab (IgG)
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.