

Product Datasheet

Estrogen Receptor beta Immunizing peptide (orb66398)

Description	Estrogen Receptor beta Immunizing peptide
Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb66398 (Estrogen Receptor beta 2 antibody)
Protein Sequence	TQDADSSRKLAH
Expiration Date	12 months from date of receipt.