


Product Datasheet

HA Tag antibody (orb66652)

Description	Mouse monoclonal to HA Tag.
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:3000-1:10000, IF:1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK 293 cell ly...