


## Product Datasheet

### S Tag antibody (orb66660)

<b>Description</b>	Mouse monoclonal to S Tag.S tag is the name of an oligopeptide derived from pancreatic ribonuclease ...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Other
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB: 1:3000-1:10000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of S-Tag fusion pr...