

Product Datasheet

IL33 antibody (orb66677)

Description	Mouse monoclonal to IL33. Interleukin 33 (IL33) is a cytokine that binds to and signals through IL1RL...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	N-terminal region of human IL33
Preservatives	0.02% Thimerosal
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Caprylic Acid-Ammonium Sulfate precipitation
MW	31
Uniprot ID	Q99JB6
Dilution Range	ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.