

Product Datasheet

PDGF Receptor beta (phospho-Tyr771) antibody (orb6668)

Description	PDGF Receptor beta (phospho-Tyr771) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PDGFRB around the phosphorylation site of Tyr771 (SN(p-Y)MA)
Target	Phospho-PDGF Receptor beta (Tyr771)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	190 kDa
Uniprot ID	P09619
Expiration Date	12 months from date of receipt.