


Product Datasheet

GADD45GIP1 antibody (orb66736)

Description	Rabbit polyclonal to GADD45GIP1. Through inhibition of cyclin-dependent kinases (CDKs), GADD45GIP1 ac...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human GADD45GIP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	25
Uniprot ID	Q8TAE8
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...