



Product Datasheet

NDUFAF4 antibody (orb667418)

Description	Rabbit polyclonal antibody to NDUFAF4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant fusion protein of human NDUFAF4. The exact sequence is proprietary.
Target	NDUFAF4
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9D1H6 , Q9NQR8 , Q9P032
Entrez	29078 , 362495 , 68493
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.




Western blot analysis of NDUFAF4 express...



Immunohistochemical analysis of NDUFAF4

...



Immunofluorescent analysis of NDUFAF4 st...