



Product Datasheet

HJV antibody (orb66746)

Description	Rabbit polyclonal to HJV.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal region of human HJV
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	Q6ZVN8
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...



IP-WB.
Immunoprecipitation (IP) was perf...