


Product Datasheet

BMP2 antibody (orb66753)

Description	Rabbit polyclonal to BMP2. Bone morphogenetic protein 2 (BMP2) is a secreted protein that belongs to ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse BMP2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	45
Uniprot ID	P21274
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...