



Product Datasheet

PHOSPHO2 antibody (orb667686)

Description	Rabbit polyclonal antibody to PHOSPHO2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human PHOSPHO2. The exact sequence is proprietary.
Target	PHOSPHO2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8TCD6 , Q66HC4 , Q9D9M5
Entrez	73373 , 295663 , 493911
Dilution Range	WB: 1-500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of PHOSPHO2 expres...



Immunofluorescent analysis of PHOSPHO2 s...