



Product Datasheet

PGC1 alpha antibody (orb667759)


Description	Rabbit polyclonal antibody to PGC1 alpha
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide of human PGC1 alpha. The exact sequence is proprietary.
Target	PPARGC1A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O70343 , Q9UBK2 , Q9QYK2
Entrez	19017 , 83516 , 10891
Dilution Range	WB: 1-500-2000, IF: 1-50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of PGC1 alpha expr...



Immunohistochemical analysis of PGC1 alp...



Immunofluorescent analysis of PGC1 alpha...