


Product Datasheet

FPN antibody (orb66800)

Description	Rabbit polyclonal to FPN.Solute carrier family 40 (iron-regulated transporter), member 1 (FPN) is a ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse FPN
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	63
Uniprot ID	Q9JHI9
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...