


Product Datasheet

MLL5 antibody (orb66810)

Description	Rabbit polyclonal to MLL5.MLL5 is a member of the myeloid/lymphoid or mixed-lineage leukemia (MLL) f...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MLL5
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	205
Uniprot ID	Q8IZID2
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...