


Product Datasheet

CD40L antibody (orb66812)

Description	Rabbit polyclonal to CD40L.CD40 ligand (CD40L) is expressed on the surface of T cells and it regulat...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal region of human CD40L
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	29
Uniprot ID	P29965
Dilution Range	WB:1:1000-1:5000,ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant hum...