


Product Datasheet

FLT3L antibody (orb66818)

Description	Rabbit polyclonal to FLT3L.Fms-related tyrosine kinase 3 ligand (FLT3L) is a cytokine that stimulate...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	C-terminal region of mouse FLT3L
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	26
Uniprot ID	P49772
Dilution Range	WB:1:1000-1:5000,ELISA:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant mouse FLT3L.