



Product Datasheet

STK4 antibody (orb668291)

Description	Rabbit polyclonal antibody to STK4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	Recombinant fusion protein of human STK4. The exact sequence is proprietary.
Target	STK4
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9JI11 , Q13043
Entrez	58231 , 6789
Dilution Range	WB: 1-500-2000, IHC-P: 1-50-200, IP: 1-20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of STK4 expression...



Immunohistochemical analysis of STK4 sta...



Immunofluorescent analysis of STK4 stain...