

Product Datasheet

Cynomolgus FOLR2 Protein (orb668875)

Description	Cynomolgus FOLR2, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 34 - His 24...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FOLR2
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.5 kDa. The protein migrates as 33-43 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Thr 34 - His 245
Purity	90%
Source	Cynomolgus FOLR2, His Tag (orb668875) is expressed from human 293 cells (HEK293). It contains AA Thr 34 - His 245 (Accession # A0A2K5U027-1).
MW	26.5 kDa
Uniprot ID	A0A2K5U027
Expiration Date	6 months from date of receipt.