
Product Datasheet

Cynomolgus IL-31 Protein (orb668883)

Description

Cynomolgus IL-31, His Tag is expressed

AA Leu 27 - Thr 16...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-31

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 17.3 kDa. The protein migrates as 20-25 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus IL-31, His Tag (orb668883) is expressed from human 293 cells (HEK293). It contains AA Leu 27 - Thr 163 (Accession # G7PJ78-1).

MW

17.3 kDa

Uniprot ID

G7PJ78

Expiration Date

6 months from date of receipt.