
Product Datasheet

Cynomolgus Mesothelin / MSLN (296-580) Protein (orb668888)

Description

Cynomolgus Mesothelin (296-580),

293 cells (HEK293). It contains AA ...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Mesothelin

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 33.7 kDa. The protein migrates as 35-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Asp 296 - Gly 580

Purity

95%

Source

Cynomolgus Mesothelin (296-580), His Tag (orb668888) is expressed from human 293 cells (HEK293). It contains AA Asp 296 - Gly 580 (Accession # UPI000C7396BF).

MW

33.7 kDa

Uniprot ID

UPI000C7396BF

Expiration Date

6 months from date of receipt.