
Product Datasheet

Mouse Mesothelin / MSLN (298-600) Protein (orb668889)

Description

Mouse Mesothelin (298-600), His Tag

(HEK293). It contains AA Asp 2...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Mesothelin

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 36.0 kDa. The protein migrates as 42-48 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_061345.1

Purity

90%

Source

Mouse Mesothelin (298-600), His Tag (orb668889) is expressed from human 293 cells (HEK293). It contains AA Asp 298 - Ser 600 (Accession # Q61468-1).

MW

36.0 kDa

Uniprot ID
[Q61468](#)
NCBI
[NP_061345.1](#)
Expiration Date

6 months from date of receipt.