

Product Datasheet

Human NRG1 Beta 1 Protein (orb668892)

Description	Human NRG1 Beta 1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ser 2 - Lys 246...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	NRG1 Beta 1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 53.2 kDa. The protein migrates as 65-80 kDa and 90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_039250.2
Purity	90%
Source	Human NRG1 Beta 1, Fc Tag (orb668892) is expressed from human 293 cells (HEK293). It contains AA Ser 2 - Lys 246 (Accession # Q02297-6).
MW	53.2 kDa
Uniprot ID	Q02297
NCBI	NP_039250.2
Expiration Date	6 months from date of receipt.