

Product Datasheet

HCoV-NL63 S1 protein (orb668902)

Description	HCoV-NL63 S1 protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Cys 19 - Va...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Spike S1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 79.5 kDa. The protein migrates as 120-140 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	Cys 19 - Val 717
Purity	90%
Source	HCoV-NL63 S1 protein, His Tag (orb668902) is expressed from human 293 cells (HEK293). It contains AA Cys 19 - Val 717 (Accession # Q6Q1S2-1).
MW	79.5 kDa
Uniprot ID	Q6Q1S2
Expiration Date	6 months from date of receipt.