


Product Datasheet

Human FGFR2b protein (orb668919)

Description	Recombinant human Fibroblast Growth Factor Receptor 2 protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the C-terminus
Note	For research use only
Protein Sequence	22-378
Purity	> 95%
Source	HEK293 cells
MW	41.2
Uniprot ID	P21802
Expiration Date	6 months from date of receipt.



Human FGFR2b protein