


## Product Datasheet

### Human FGFR3c protein (orb668922)

<b>Description</b>	Recombinant human Fibroblast Growth Factor Receptor 3 protein
<b>Reactivity</b>	Human
<b>Tested Applications</b>	Functional Studies, IA, WB
<b>Immunogen</b>	Unconjugated
<b>Preservatives</b>	PBS 20% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20° or -70°C
<b>Tag</b>	10xHis at the C-terminus
<b>Note</b>	For research use only
<b>Protein Sequence</b>	23-375
<b>Purity</b>	> 95%
<b>Source</b>	HEK293 cells
<b>MW</b>	39.7
<b>Uniprot ID</b>	<b>P22607</b>
<b>Expiration Date</b>	6 months from date of receipt.



Human  
FGFR3c  
protein