


Product Datasheet

Human Ecto-5'-nucleotidase, 5'-NT or CD73 protein (orb668931)

Description	Recombinant human Ecto-5'-nucleotidase, 5'-NT or CD73 protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	6xHis at the C-terminus
Note	For research use only
Protein Sequence	27-549
Purity	> 95%
Source	HEK293 cells
MW	58.8
Uniprot ID	P21589
Expiration Date	6 months from date of receipt.



Human Ecto-5'-nucleotidase, 5'-NT or CD7...