



Product Datasheet

TICAM1 antibody (orb668962)

Description	Goat polyclonal antibody to TICAM1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Target	TICAM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Small Intestine. Recommended concentration: 2.5µg/ml.. Notes: Flow Cytometry: Flow cytometric analysis of K562 cells. Recommended concentration: 10ug/ml.
Protein Sequence	HARADEHIALRVREK
Clonality	Polyclonal
MW	76.4
Entrez	148022
Dilution Range	ELISA: 1:64000, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Sm...



Flow cytometric analysis of paraformalde...