
Product Datasheet

Cynomolgus IL-1 Rrp2 / IL-1 R6 Protein (orb669283)

Description

Cynomolgus IL-1 Rrp2, His Tag is

(HEK293). It contains AA Asp 20 - Ar...

Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-1 Rrp2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.9 kDa. The protein migrates as 50-66 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Purity	90%
Source	Cynomolgus IL-1 Rrp2, His Tag (orb669283) is expressed from human 293 cells (HEK293). It contains AA Asp 20 - Arg 335 (Accession # EHH61702.1).
MW	37.9 kDa
NCBI	EHH61702.1
Expiration Date	6 months from date of receipt.