

Product Datasheet

Cynomolgus CD24 Protein (orb669289)

Description	Cynomolgus CD24, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr 27 - Lys 54 (...)
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD24
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 29.3 kDa. The protein migrates as 42-50 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	Thr 27 - Lys 54
Purity	90%
Source	Cynomolgus CD24, Fc Tag (orb669289) is expressed from human 293 cells (HEK293). It contains AA Thr 27 - Lys 54 (Accession # I7GKK1-1).
MW	29.3 kDa
Uniprot ID	I7GKK1
Expiration Date	6 months from date of receipt.