

Product Datasheet

Human LAIR1 protein (orb669302)

Description	Recombinant Human Leukocyte-Associated Immunoglobulin-Like Receptor 1 is produced by human 293 cells...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 2umM PB, 15umM NaCl, pH 7.2
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	Q6GTX8
Expiration Date	6 months from date of receipt.