


Product Datasheet

PRMT10 antibody (orb66931)

Description	Rabbit polyclonal to PRMT10. Protein arginine methyltransferase 10 (PRMT10) belongs to the protein ar...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human PRMT10
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	95
Uniprot ID	Q6P2P2
Dilution Range	WB:1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extr...