

Product Datasheet

Mouse B7-2/CD86 protein (orb669311)

Description	Recombinant Mouse T-lymphocyte Activation Antigen CD86 is produced by our Mammalian expression syste...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Uniprot ID	P42082
Expiration Date	6 months from date of receipt.