

## Product Datasheet

### Human B7-2/CD86 protein (orb669313)

<b>Description</b>	Recombinant Human CD86 is produced by our Mammalian expression system and the target gene encoding A...
<b>Endotoxins</b>	Less than 0.1 EU/ug
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P42081</b>
<b>Expiration Date</b>	6 months from date of receipt.