

Product Datasheet

SARS-CoV-2 Spike (RBD) protein (orb669319)

Description	Recombinant SARS-CoV-2 Spike Protein is produced by HEK293. The target gene encoding 319-532 is expr...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-hFc
Note	For research use only
Uniprot ID	PODTC2
Expiration Date	6 months from date of receipt.