

Product Datasheet

Mouse Interferon alpha-1 protein (orb669396)

Description	Recombinant Mouse Interferon alpha-1 is produced by Human 293 Cells. The target gene encoding C24-K1...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-8His
Note	For research use only
Uniprot ID	P01572
Expiration Date	6 months from date of receipt.