

Product Datasheet

Mouse PDGF-B protein (orb669401)

Description	Recombinant Mouse PDGF-B is produced by E.coli. The target gene encoding S82-T190 is expressed with ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, um Urea pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	N-6His
Note	For research use only
Uniprot ID	P31240
Expiration Date	6 months from date of receipt.