

Product Datasheet

Human ACE-2 protein (orb669405)

Description	Recombinant Human Angiotensin-Converting Enzyme 2 is produced by our Mammalian expression system and...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 300mM NaCl, 10% Glycerol, 100mM Glycine, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-Fc
Note	For research use only
Uniprot ID	Q9BYF1
Expiration Date	6 months from date of receipt.