

Product Datasheet

Human b-Lactamase-Like 2 protein (orb669422)

Description	Recombinant Human beta-Lactamase-Like Protein 2 is produced by our E.coli expression system and the ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS,1mMDTT pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-GST
Note	For research use only
Uniprot ID	Q53H82
Expiration Date	6 months from date of receipt.