

Product Datasheet

Human Serine Protease HTRA2 protein (orb669497)

Description	Recombinant Human High Temperature Requirement Protein-2 is produced by our Mammalian expression sys...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	O43464
Expiration Date	6 months from date of receipt.