

Product Datasheet

Human Carbonic Anhydrase-Related 11 protein (orb669501)

Description	Recombinant Human Carbonic Anhydrase 11 is produced by our Mammalian expression system and the target...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	O75493
Expiration Date	6 months from date of receipt.