

Product Datasheet

Human B3GNT1 protein (orb669508)

Description	Recombinant Human N-Acetylglucosaminidase beta-1,3-N-Acetylglucosaminyltransferase is produced by our ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris HCl,150mM NaCl,Ph7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Uniprot ID	O43505
Expiration Date	6 months from date of receipt.