

Product Datasheet

Human Peptidyl-Prolyl Cis-Trans Isomerase D protein (orb669547)

Description	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase D is produced by our E.coli expression system ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 10% glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&C-His
Note	For research use only
Uniprot ID	Q08752
Expiration Date	6 months from date of receipt.