

## Product Datasheet

# Human Ubiquitin-Conjugating Enzyme Variant 2 protein (orb669561)

<b>Description</b>	Recombinant Human Ubiquitin-Conjugating Enzyme E2 Variant 2 is produced by our E.coli expression sys...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 50mm HEPES, 150mM NaCl, pH 7.0.
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	N-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q15819</a>
<b>Expiration Date</b>	6 months from date of receipt.